

## Voice Selector Converse and wireless earphones

### Safety Instructions

#### How to use the device:

Voice Selector Converse including a wireless earphones is a unique assistive listening device, reducing ambient noise and focusing audio in the selected direction.

Storage and operation temperatures should be between 10° and 45°C, and humidity 10% to 95%.

#### How to clean the device:

The device should be cleaned with a damp cloth.

#### Instructions for disposing of electronic equipment:

Electronic waste must be treated according to local regulations. Nuance Hearing Ltd. recommends recycling the device at the regional disposal of electronic equipment. Local laws must be followed.

#### According to the regulations, we are required to explain the following warnings:

- Warning: Prolonged listening at high volume is not recommended and may cause irreversible damage to your hearing
- Warning: Use of accessories, components or spare parts not supplied by Nuance Hearing may impact use and void the warranty
- The battery in Voice Selector Converse and in the wireless earphones are not replaceable and should not be replaced

#### Standards and regulations

Voice Selector Converse and the wireless earphones meets the strictest international standards for electrical safety and electromagnetic compatibility.

